


# Advantage Pole 350™ / Portable Options


A solid, floor mounted pole suitable for users who require moderate assistance to achieve small steps during transfers.

## Features and Benefits


- **Flexible** - The Advantage Pole 350™ is ideal for applications where the ceiling height / structure will not accommodate a SuperPole™ (a floor to ceiling grab bar). The Advantage Pole 350™ mount directly to the floor, so they can be installed virtually anywhere.
- **Portable Options** - Each portable model quickly secures and releases from a low profile floor mounted plate. Additional floor plates can be purchased to use the Advantage Pole™ in multiple locations.



Bedroom



Bathroom



Hospital



Advantage Pole™ Portable

## Technical Specifications:

Advantage Pole™	54" / 137 cm high.
Pole Mast	2.25" / 5.7 cm diameter and 28" / 71 cm high.
Floor Plate	Six hole floor screw pattern (6 X 3/8" holes on 6.75" diameter)
Pole Material	Heavy gauge steel, 1.5" / 3.81 cm diameter.
	High temperature baked powder coat gloss white enamel.
Grip Material	Grey, premium quality, 20" / 51 cm long.
	Does not contain natural dried rubber or latex.
Portable Floor Plate	Six hole floor screw pattern (6 X 3/8" holes on 6.75" diameter)
Weight Capacity	350 lbs / 159 kg.
Assistance Level	

