


# Advantage Rail™ / Portable Options


A solid, floor mounted pivoting and locking bar suitable for users who require moderate assistance to achieve small steps during transfers.

## Features and Benefits


- **Maximum Support** - The Advantage Rail™ has a pivoting feature that allows the horizontal rail to be positioned where it can provide maximum support. The rail locks in an instant in case there is a loss of balance, a need for rest, or a desire to get to a seated position.
- **Flexible** - The Advantage Rail™ is ideal for applications where the ceiling height / structure will not accommodate a SuperPole™ or SuperPole™ with SuperBar™ (a floor to ceiling grab bar). The Advantage Rail™ and Advantage Rail 350™ mount directly to the floor, so they can be installed virtually anywhere.
- **Portable Options** - Each portable model quickly secures and releases from a low profile floor mounted plate. Additional floor plates can be purchased to use the Advantage Rail™ in multiple locations.
- Aluminum Advantage Rail "T" Handle available. Weight capacity 160 lbs.



Bedroom



Bathroom



Hospital



Advantage Rail™ Portable

## Technical Specifications:

Advantage Rail™	Adjustable range 30"-38" / 75-95 cm. 18" / 48 cm long with premium quality grey grip.
Rail Mast	2.25" / 5.7 cm diameter and 28" / 71 cm high.
Floor Plate	Six hole floor screw pattern (6 X 3/8" holes on 6.75" diameter)
Rail Material	Heavy gauge steel, 1.5" / 3.81 cm diameter.
	High temperature baked powder coat gloss white enamel.
Grip Material	Grey, premium quality, 20" / 51 cm long.
	Does not contain natural dried rubber or latex.
Portable Floor Plate	Six hole floor screw pattern (6 X 3/8" holes on 6.75" diameter)
Weight Capacity	250lbs / 113 kg. Advantage Rail 350™ weight capacity 350 lbs / 159 kg.
Assistance Level	

