


SuperPole™ with Safety Island™

A unique circular ring with a tray that mounts at any height to the SuperPole™. Combines the fixed support benefits of a grab bar with the simplicity and convenience of a tray.

Features and Benefits

- **360 Degree Support** – Effective ring design provides continuous support to make small, safe steps, reducing unnecessary reaching and twisting. Ring & tray will not spin freely around the pole, providing fixed support and an ideal weight bearing surface when gripping strength is weak.
- **Convenient** – Easy to clean tray can be accessed from around the SuperPole™. Perfect for reading material, medication or items needed within arms reach. The Safety Island™ will mount to the SuperPole™ at virtually any height.
- **Ideal for Kitchens** – Install in the center of the kitchen to create secure support when moving between chair, table, sink or counter. Use the tray to place meals while moving from cooking surface to kitchen table.



Technical Specifications:

Tubular Ring	Typically 28-38"/71-97cm from floor, 18"/46cm diameter, 1.25"/3.175cm diameter tube for easy gripping.
Tray	14"/36cm diameter A.B.S. tray, easy to wipe clean. Weight capacity 15lbs/7kg.
Floor to Ceiling Range	93-99"/236-251cm (standard).
Uni-Fit Extender™	Increase floor to ceiling range to 100-120"/254-305cm. Sold separately. NOT TO BE USED WITH SUPERTRAPEZE™.
SuperPole™ Ultra	NOT TO BE USED WITH SAFETY ISLAND™.
Ceiling Plate Extender™	Allows standard ceiling plate (20"/51cm length) to span 30"/76cm with holes for ceiling attachment to 24"/61cm ceiling beam centers. Sold separately.
Pole Material	Heavy gauge steel, 1.5"/3.81cm diameter. High temperature baked powder coat gloss white enamel.
Grip Material	Grey, premium high quality grip, 20"/51 cm long. Does not contain natural dried rubber or latex. Washable.
Expansion Mechanism	Steel jackscrew, zinc plated, 1-1/8" fine cut thread.
Weight Capacity	250lbs/113kg. Bariatric model 450lbs/205kg capacity.
Assistance Level	

