


# SuperPole™ with SuperTray™

The SuperTray™ easily snaps onto the SuperPole™ to provide the convenience of a bedside or overbed table, as well as transfer support from the pole.

## Features and Benefits

- **Flexible** - The tray is height adjustable and easily installed at the most desirable and convenient location for the user.
- **Tray** - The tray has a generous flat surface perfect for meals, medications, reading material or anything else one may need within arms reach. A convenient surface for all bedside needs.
- **Safety** - Tray surface comes complete with non-skid mat to prevent items from slipping. A nice precaution when handling hot liquids.
- **Versatile** - The SuperTray™ easily pivots out of the way when not in use, it's slim profile allows for easy storage and cleaning.



## Technical Specifications:

SuperTray™	Generous 11" x 15" (28 x 38 cm) A.B.S.
	SuperTray™ and SuperBar™ cannot be mounted at the same height.
Floor to Ceiling Range	93" to 99" / 236 to 253 cm (standard).
Uni-Fit Extender™	Increase floor to ceiling range from 100" to 120" / 254 to 305 cm.
SuperPole™ Ultra	Accommodates floor to ceiling range from 121"-140" / 307 to 356 cm.
Ceiling Plate Extender™	Allows standard ceiling plate 20" / 51 cm length to span 30" / 76 cm with holes for ceiling attachment at 24" / 61 cm ceiling beam centers.
Pole Material	Heavy gauge steel, 1.5" / 3.81 cm diameter. High temperature baked powder coat gloss white enamel.
Grip Material	Grey, premium quality, 20" / 51 cm long. Does not contain natural dried rubber or latex.
Expansion Mechanism	Steel jackscrew, zinc plated, 1-1/8" fine cut thread.
Weight Capacity	250lbs / 113 kg for SuperPole™. Plastic tray holds up to 15 lbs / 7 kg.
Assistance Level	



NOTE: The SuperTray™ is NOT to be used as a transfer aid. Assistance Level pertains to the SuperPole™.

