

SPS

Single Patient Slings



UD SERIES

Apex *Single Patient Slings* have been specially designed to meet the needs of today's healthcare environments. These cost effective slings are patient specific and will greatly reduce the incidence of cross contamination found in acute care and other facilities . The Apex *Single Patient Slings* is made from a comfortable yet durable non-woven material. Control is achieved by placing the name of the patient in the designated box on each label. Available in a variety of sizes and styles the Apex *Single Patient Slings* may be safely used with most patient lift systems. Reduce costs and increase patient health and safety with Apex SPS today!



UDH SERIES

- Lower cost alternative to regular slings
- For use with one patient only
- 600lb safe working load
- Do not wash (discard immediately if washed)
- Packed 10 per carton and individually poly wrapped
- Color coded for size identification
- Recommended for use by professional caregivers

SPECIFICATIONS:



SL-UD812
SL-UD813

Universal split leg medium
Universal split leg large

SL-UDH821
SL-UDH822
SL-UDH823

Universal small with head support
Universal medium with head support
Universal large with head support

