

Spryte STAND AID LIFT



The Spryte stand aid is a transport assist unit which keeps the resident active and engaged in the process. Transport functions of all types are quick and require minimal caregiver assistance. Users grasp the middle bar and using their own strength pull themselves up into position. A padded split seat swings out to allow loading or unloading then swings back to form a comfortable and secure seat for transport. An excellent alternative to a wheelchair for easier toileting. Note: Spryte is designed for residents that have the strength and stability to lift and support themselves. For residents with limited strength and stability an electric stand assist is recommended.

STA182



- 400 lb safe working load
- Dual knee pads for secure transports
- Sturdy yet easily fits through narrow doorways
- Affordable alternative to electric stand aids



Optional Safety Strap



Spryte Stand Aid Lift Specifications

Maximum width	Maximum height	Base height	Overall length	Unit weight
25.5"	43"	4.5"	32"	60 lbs.

ISO
9001
CE